

## Vitrification Straw

**Cryoteb** is the special vitrification container consisting of a thin film strip attached to a hard plastic handle.

### Main Features

- 10 different colors.
- Allows up to 4 specimens to be loaded per device
- High transparency.
- Black mark to know on which side the specimens must be loaded.
- Manufactured in malleable plastic resistant to liquid nitrogen.
- Very cost effective

### QUALITY CONTROL

- . MEA (Mouse Embryo Assay): One cell assay  $\geq 80\%$  after 96 hours
- . Endotoxin:  $< 20\text{EU/device}$
- . Sterility Test (Bacteria, Fungi)
- . Cytotoxicity Test
- . Sensitization Test
- . HSSA (Human sperm survival assay):  $\geq 70\%$  motility at 24 hours



2.5 Euro