


ICSI PIPETTE

For intracytoplasmic injection of sperm inside the oocyte.



- It is easy to control the sperm with long arm.
- Short spike minimize damage to cell cytoplasm.
- Appearance of respective pipette are inspected one by one with image system.

Quality control test:

- Mouse embryo assay (MEA) ≥ 85
- Endotoxin ≤ 0.5 EU/device
- SAL (Sterility Assurance Level)
- Irritation test – non-irritate
- Cytotoxicity – non cytotoxic

Shelf life:

- 3 years

Storage condition:

- Dry and cool place (15-25 °)

	ID (μM)	OD (μM)	Angle
ICSI	4.5	6	30