

## Vitrification straw (CRYOTEB)

### Short handle model

Cryoteb is the Embryo / Oocyte vitrification device

Cryoteb is used for storage of vitrified oocytes and embryos.



- Vitrification straw device is consisting of a fine strip of transparent film attached to a plastic handle resistant to liquid nitrogen.
- Available in 10 different colors.
- Cryoteb allows up to 4 specimens to be loaded per device.
- Cap covering to protect the transparent film.
- Cap with steel prevent floating in liquid nitrogen.

#### Quality control test:

- Mouse embryo assay (MEA)  $\geq 85$
- Endotoxin  $\leq 0.5$  EU/device
- SAL (Sterility Assurance Level)
- Irritation test – non-irritate
- Cytotoxicity – non cytotoxic

#### Shelf life:

- 3 years
- Storage condition:
- Dry and cool place (15-25 °)

