

EMBRYO TRANSFER CATHETER

Used to place in vitro fertilized (IVF) embryos into the uterine. Two model of ET catheter is available: simple & Echo tip model.



- ✓ The guide catheter is curved to facilitate insertion.
- ✓ The guide catheter's rounded bulb tip eases passage through the cervix and reduce the trauma during the insertion.
- ✓ Their 30° curved guide and characteristic flexibility allow easy navigation through the channel of the cervix in any kind of anatomy.
- ✓ The transfer catheter is so soft and flexible.
- ✓ The special design minimizes the volume of media that is required.
- ✓ The cervical stop is set at 4 cm and can be adjusted to 5 cm.
- ✓ The catheters are packaged in two packaging, hard vacuum and sterile pouch.

Quality control test:

Shelf life:

- 3 years
- Storage condition:
- Dry and cool place (15-25 °)
- Mouse embryo assay (MEA) ≥ 85
- Endotoxin ≤ 0.5 EU/device
- SAL (Sterility Assurance Level)
- Irritation test – non-irritate
- Cytotoxicity – non cytotoxic

