

BIOPSY PIPETTE

For removing a blastomeric cell or polar body for pre-implantation genetic diagnosis.



Quality control test:

- Mouse embryo assay (MEA) ≥ 85
- Endotoxin ≤ 0.5 EU/device
- SAL (Sterility Assurance Level)
- Irritation test – non-irritate
- Cytotoxicity – non cytotoxic

Shelf life:

- 3 years

Storage condition:

- Dry and cool place (15-25 °)

| | ID (μM) | OD (μM) | Angle |
|----------------|----------------------|----------------------|-------|
| Biopsy pipette | 15 | 20 | 30 |
| | 20 | 25 | 30 |
| | 25 | 35 | 30 |