

HOLDING PIPETTE

Holding micropipettes have been designed to provide a high level of specimen control of the oocyte, embryo or blastocyst during intracytoplasmic sperm injection or PGT techniques.



- Manufactured with special process in order to hold the oocytes flat and horizontal.

Quality control test:

- Mouse embryo assay (MEA) ≥ 85
- Endotoxin ≤ 0.5 EU/device
- SAL (Sterility Assurance Level)
- Irritation test – non-irritate
- Cytotoxicity – non cytotoxic

Shelf life:

- 3 years

Storage condition:

- Dry and cool place (15-25 °)

	ID (μM)	OD (μM)	Angle
Holding pipette	15-20	90	30
	15-20	100	30
	15-20	85	30